

المادة: جزيئي عملي / 5

المرحلة: الرابعة

اسم الطالب:

السعر (500)

## Plasmid DNA extraction

A plasmid is a super coiled (covalently closed circular CCC) extra-chromosomal DNA molecule, it is capable of replicating independently from chromosomal DNA.

- ❖ There are many methods are used fro plasmid DNA isolation, such as:
  - a) Alkaline lysis method
  - b) Boiling lysis method
  - c) Phenol method
  - d) Ethidium bromide- Caesium chloride density gradient centrifugation method
  - e) Mini- column

### Alkaline lysis protocol for plasmid extraction:

- **Harvesting**

- 1) Inoculate *E. coli* (RRI) carried PBR 322 plasmid into 2 ml of LB medium containing the appropriate **ampicillin 100 mg/ml** in 15 ml tube volume. Incubate the culture **overnight at 37C in shaking incubator.**
- 2) Pour 1.5 ml of the culture into a microfuge tube. **Centrifuge for 30 sec. at 12,000 rpm in a microfuge at 4 C.**
- 3) Remove the medium, leaving the bacterial pellet as dry as possible.

- **Lysis by alkali:**

- 1) Resuspend the bacterial pellet (obtained from step 3 above) in **100 µl of ice- cold solution I** by vigorous avertexing.

Solution I
50 mM glucose
25 mM Tris -Hcl
10 mM EDTA

Solution I can be prepared in batches of approximately 100 ml, adjust the pH to 8 with HCl then autoclave for 15 min. at 10 lb/in 2 hr. and store at 4C.


- 2) Add 200 µl of freshly prepared **solution II**, close the tube tightly and mix the contents by inverting the tube gently (**Do not vortex**). and store the tube on ice.

Solution II
0.2 N NaOH (freshly diluted from a 1 N stock)
1% SDS

- 3) Add 150  $\mu$ l of **ice – cold solution III**. Close the tube and mix gently by inversion (6 times), then disperse solution III through the viscous bacterial lysate and store the tube on ice for 3-5 min.

Solution III	
5 M potassium acetate	60 ml
Glycyl acetic acid	11.5 ml
H <sub>2</sub> O	28.5 ml

- 4) Centrifuge milky suspension (emulsion) at 12000 rpm for 5 min. in a microfuge Eppendorf tube. Transfer **the supernatant** to a new tube.
- 5) Add **450  $\mu$ l of phenol : chloroform (ratio 25/25)**. Mix by vortexing, after centrifugation at 12000 rpm for 2 min. in a microcentrifuge, transfer the aqueous phase to a new tube.



- 6) Add **(0.7 volume) of ice- cold isopropanol 100%** to supernatant, mix by vortexing, centrifuge at 12000 rpm for 20 min. at 4 C.
- 7) Rinse the pellet of double strand DNA (super coiled) with **1 ml of 70% ethanol**. Remove the supernatant and allow the pellet of nucleic acid to dry in the air for 10 min.
- 8) Re-suspend the nucleic acids in **50  $\mu$ l of TE (pH 8)** containing DNase- free RNase solution (20  $\mu$ g/ml). vortex briefly and store the DNA at 4 C.
- 9) Test the size of plasmid and finding protocol using gel electrophoresis technique.

## Extraction of RNA

Main types of RNA are:

- 1) **mRNA:** transfer the genetic information from DNA to protein synthesis place (ribosome) in cytoplasm.
- 2) **tRNA:** found in cytoplasm, carrying amino acids for arranging them on the strand of mRNA. There are more than 20 kinds of tRNA in each cell.
- 3) **rRNA:** it calls ribosomal RNA, because it builds the ribosomal structure with nucleoproteins, have a role in protein synthesis in the cytoplasm.

### To isolate RNA molecules, we have to follow the steps below:

#### A) Lysis the cell

- ❖ Treat the suspension of *E. coli* (0.5 O.D at 600 nm) with equal volume of **Diethyl pyrocarbonate (DEPE)**, which work as **inhibiter to RNase (exogenous RNase)**. Homogenized and incubate at 4 C.
- ❖ Add (0.4 ml) of 5% SDS to the cell suspension  
(SDS acts on lipid of membrane that allowed access of lysozyme to peptidoglycan, also SDS inhibits RNase).
- ❖ Use (2ml) of lysozyme (400 µg/ml) containing buffer per (10 ml) of *E. coli*, homogenized and incubated at room temperature (37C) for 5-20 min.

#### B) Purification of RNA

- ❖ Centrifuge the suspension of lysis cells at 3000 rpm for 10 min.
- ❖ Suspend the pellet with equal volume of **guanidinium thiocyanate** which acts as **strong inhibitors of RNase and denaturation of protein**. The pH should not be alkaline due to liability of rRNA.
- ❖ DNA can be removed by using **DNase I** to minimize genomic DNA contamination.

- ❖ The contaminant proteins are removed by adding a mixture of **phenol: chloroform: isoamyl alcohol (25:24:1)**, and repeat this step many times.
- ❖ Transfer the upper layer to another tube and add **1 volume of ice cold absolute isopropanol** to precipitate RNA.
- ❖ Centrifuge 1000 rpm for 20 min. at 4 C
- ❖ Discard supernatant and wash pellet with **70% ethanol**
- ❖ Centrifuge 1000 rpm for 5 min. at 4 C
- ❖ Re-suspend the RNA in **100 µl of TE buffer** and frozen

### Calculate the Concentration of RNA:

$$\text{Concentration of RNA } (\mu\text{g/ml}) = \text{O.D } 260 \text{ nm} \times 40 \times \text{dilution factor}$$

- 1 O.D 260 nm = 40 µg/ml of RNA

### Calculate the Purity of RNA:

$$\text{O.D } 260 \text{ nm} / \text{O.D } 280 \text{ nm} = \text{should be } 2$$